PROTANK-1 323-A

Ca²⁺ 1.5 mmol/L (3 mEq/L)

K⁺ 2 mmol/L (2 mEq/L)

Glucose 1 g/L

EN Acid concentrate solution for hemodialysis and related therapies
FR Solution acide concentrée pour hémodialyse et les thérapies connexes
ES Solución de concentrado ácido para hemodiálisis y terapias relacionadas
PT Solução ácida concentrada para hemodiálise e terapias relacionadas
IT Soluzione acida concentrata per emodialisi e terapie correlate

NL Zuurconcentraat oplossing voor hemodialyse en aanverwante therapieën
BS Koncentrat kiseline za hemodijalizu i srodne terapije
HR Kiseli koncentrat za hemodijalizu i srodne terapije
DE Säurekonzentratlösung für die Hämodialyse und verwandte Therapien

Composition • Composition • Composición • Composição • Composizione • Samenstelling • Sastav • Sastav • Zusammensetzung

NaCl: 270.88 g/L

MgCl₂ 6 H₂O: 4.57 g/L

CH₃COOH: 8.11 g/L

H₂O Ph. Eur. q.s.: 1000 mL

CaCl₂ 2 H₂O: 9.92 g/L

KCl: 6.71 g/L

C₆H₁₂O₆ H₂O: 49.50 g/L

EN Electrolytic composition of the ready-to-use solution
FR Composition électrolytique de la solution prête à l'emploi
ES Composición electrolítica de la solución lista para su uso
PT Composição eletrolítica da solução pronta a usar
IT Composizione elettrolitica della soluzione pronta all'uso
NL Elektrolytische samenstelling van de gebruiksklare oplossing

BS Elektrolitni sastav rastvora spremnog za upotrebu
HR Sastav elektrolita u otopini spremnoj za korištenje
DE Elektrolytische Zusammensetzung der gebrauchsfertigen Lösung

mmol/L:

1 L (323-A): 1.4 L (B): 42.6 L (H:O)

Na ⁺	Ca ²⁺	Mg ²⁺	K ⁺	Cl ⁻	CH ₃ COO ⁻	HCO ₃ ⁻	C ₆ H ₁₂ O ₆	Osmolarity
140	1.5	0.5	2	109	3	34	5.55	295.5 mOsm/L

Important: Store at a temperature no lower than 4°C. Pyrogen-free product.
Intended purpose: Preparation of dialysis fluid used in hemodialysis and related therapies. **Indication:** Treatment of patients with acute or chronic renal failure when hemodialysis is prescribed by a physician. **Clinical Benefits:** Facilitates the clearance of blood from uremic toxins that are able to pass through the dialysis membrane and balances out the electrolyte imbalance. **Contraindications:** There are no special contraindications for the use of NiproTANK-1. Patients who are hypersensitive to any components used in the NiproTANK-1 should not be dialyzed with this product. **Side effects:** There are no side effects reported for NiproTANK-1. **Instructions:** Acid concentrate with dilution 1:44 to be used in combination with powder bicarbonate concentrate and hemodialysis machines. To be connected to the acid concentrate connector of the central concentrate supply system. The acid concentrate from the central concentrate supply system must be used with a dialysis machine, following the instructions of the machine manufacturer. Dilution is based on the following proportion: See the product label. B is the bicarbonate concentrate and H2O is deionized water, complying with EP, ISO 23500-5, and applicable local quality requirements. **Warnings:** To be used by healthcare professionals. Home dialysis is excluded. To be used under the physician's order. Ensure that the composition corresponds to the physician's prescription. Do not use if the package is damaged. The concentrated solution must be diluted immediately before use and opened, the product must be used continuously without disconnecting from the central concentrate supply system. Do not use after expiry date stated on the package. See table below for an indication of expected dialysis fluid volume obtained from each container volume. Any unusual incident that has occurred in relation to the product should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is based. Dispose of the product according to local and national regulations.

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Theramed AG Therapie und Medizintechnik
UK Responsible Person: Nipro Medical UK Ltd.