

ACID CONCENTRATES

BAGS, CANISTERS, OCTABINS



ACID CONCENTRATES

BAGS, CANISTERS, OCTABINS

Ready-to-use liquid acid concentrates intended for use with bicarbonate concentrate.

Acid concentrate is an essential component in the preparation of dialysis fluid, whereby the hemodialysis machine mixes liquid acid concentrate with bicarbonate concentrate and purified water.

The standard mixing ratio of Nipro acid concentrates is 1:44.



PORTFOLIO

	Na ⁺	Ca ⁺⁺	Mg ⁺⁺	K ⁺	Cl ⁻	CH ₃ COO ⁻	HCO ₃ ⁻	Glucose	Osmolarity
	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	g/l	mOsm/l
306-A (A+B)	139	1,25	0,48	1,50	103,14	2,85	37,96	1,50	294,50
307-A (A+B)	139	1,75	0,48	1,50	104,14	2,85	37,96	1,50	296,00
308-A (A+B)	139	1,50	0,48	1,50	103,64	2,85	37,96	1,00	292,48
310-A (A+B)	139	1,50	0,50	1,00	103,20	2,85	37,96	1,00	291,56
313-A (A+B)	140	1,50	0,50	1,50	108,50	3,00	34,00	1,00	294,55
314-A (A+B)	140	1,25	0,50	1,50	108,00	3,00	34,00	1,00	293,80
316-A (A+B)	140	0,63	0,25	1,50	102,44	4,00	36,81	1,00	291,18
318-A (A+B)	140	1,75	0,50	2,00	109,50	3,00	34,00	1,50	299,08
319-A (A+B)	139	1,75	0,48	3,00	105,64	2,85	37,96	1,50	299,00
320-A (A+B)	139	1,25	0,50	3,00	104,70	2,85	37,96	1,50	297,59
321-A (A+B)	140	1,75	0,50	2,00	109,50	3,00	34,00	1,00	296,30
323-A (A+B)	140	1,50	0,50	2,00	109,00	3,00	34,00	1,00	295,55
324-A (A+B)	140	1,25	0,50	2,00	108,50	3,00	34,00	1,00	294,80
326-A (A+B)	139	1,50	0,48	2,50	104,64	2,85	37,96	1,50	297,25
327-A (A+B)	140	1,50	0,50	3,00	110,00	3,00	34,00	1,00	297,55
330-A (A+B)	140	1,75	0,50	3,00	110,50	3,00	34,00	1,00	298,30
331-A (A+B)	140	1,25	0,50	3,00	109,50	3,00	34,00	1,00	296,80
332-A (A+B)	140	1,50	0,50	4,00	111,00	3,00	34,00	1,00	299,55

PACKAGING DETAILS

	Volume (L)	Reference code	Boxes per pallet	Units per box	Pallet layers	Units per pallet
Bags - Regular box	3,8	BSOL45-***A4L3BX	48	3	6	144
	4,7	BSOL45-***A5LBXEU	64	2	8	128
Bags - Big box	3,8	BSOL45-***A4LBBEU	4	40	2	160
	4,7	BSOL45-***A5LBBEU	4	35	2	140
Canisters - Stackable	3,9	SOL45-***A4L	N/A	N/A	4	132
	5	SOL45-***A5L	N/A	N/A	4	132
	10,5	SOL45-***A11L	N/A	N/A	3	60
Canisters - Regular box	3,9	SOL45-***A4LBX	64	2	4	128
	5	SOL45-***A5LBX	64	2	4	128
Canisters - Big box	3,9	SOL45-***A4LBB	1	84	1	84
	5	SOL45-***A5LBB	1	84	1	84
	10,5	SOL45-***A11LBB	1	40	1	40
Octabin - Bottom connector	500	CSOL45-***A500L	N/A	N/A	1	1
Octabin - Top connector	300	CSOL45-***A300LT	N/A	N/A	1	1
	500	CSOL45-***A500LT	N/A	N/A	1	1
	600	CSOL45-***A600LT	N/A	N/A	1	1

* acid concentrate formula

Not all products are registered and approved for sale in all countries or regions. Indications of use may also vary by country and region. Please contact your Nipro representative for product availability and information.

SHELF LIFE: Bags and Canisters: 3 years | Octabin: 2 years

Nipro Renal Care is part of Nipro Corporation Japan, a leading global healthcare company established in 1954. With over 35.000 employees worldwide, Nipro serves the Medical Device, Pharmaceutical, and Pharmaceutical Packaging industries.

Nipro Renal Care is a global market leader with over 6 decades providing renal solutions for dialysis and dialysis-related treatment. We specialize in developing dialysis machines, water treatment systems, and a comprehensive portfolio of disposable medical equipment.

In order to address the needs of patients, healthcare professionals, and procurement managers alike, Nipro Renal Care is driven by innovation and patient safety to offer the highest quality products that optimize time, effort, and costs.

BECAUSE EVERY LIFE DESERVES AFFORDABLE CARE



www.nipro-group.com/en-en/our-company/our-locations

Please contact your local representative for more information.

